

# T4 DNA Ligase

**Catalog:** E-102

**Lot:** 3

**Quantity Supplied:** 10 mg in 552 ul storage buffer.

**Concentration:** 18.1 ug/ul (highly concentrated)

**Storage Buffer:** 100 mM Tris-HCl, pH 8.0, 100 mM NaCl, 1.0 mM EDTA, 200 mM imidazole, 10 mM dithiothreitol, 0.10% Triton X-100, 25% glycerol.

**Purity:** Greater than 98% as determined by SDS-PAGE gel analysis.

**DNA Ligase Activity:** approximately 600 cohesive end ligation units per ug protein using a qualitative blunt end dsDNA ligase assay (the enzyme activity per ug is roughly the same as the activity of New England Biolabs T4 DNA ligase, which was used for comparison).

**Storage:** -20 °C or below. Repeated freeze-thaw cycles have no effect on this product.

**Description:** T4 DNA ligase ligates adjacent 3'-hydroxyl and 5'-phosphate bases in double stranded DNA. T4 DNA ligase also ligates blunt ended dsDNA fragments. The recommended 10X buffer is 500 mM Tris-HCl, pH 7.5, 50 mM MgCl<sub>2</sub>, 50 mM dithiothreitol, 10 mM ATP.

**Recommended Use:** This is a highly concentrated preparation of T4 DNA ligase. Typical ligation reactions can be easily accomplished by using a 100 fold (or more) dilution of this T4 DNA ligase preparation. The reaction should be extracted with an equal volume of phenol:CHCl<sub>3</sub>, prior to agarose gel electrophoresis, in order to avoid electrophoretic artifacts.

**Enzyme Stability:** This enzyme preparation should be stored in a freezer, at -20 °C or lower in order to preserve enzyme activity long term. Stability experiments show that storage of this stock solution at 25 °C (room temperature) for 6 days results in a loss of about 25% of ligase activity. Storage at 37 °C for 6 days results in a loss of about 60% of ligase activity.

## Quality Control Tests

(1) Endonuclease contamination: None detected.

Incubation of 54 ug (~ 32,000 units) with 15 ug of pUC19 supercoiled plasmid in 50 ul 50 mM Tris-HCl, pH 7.5, 5 mM MgCl<sub>2</sub>, 5 mM dithiothreitol, at 37 °C for 3 hours resulted in no observable conversion of supercoiled plasmid to open circular form.

(2) Single Stranded Exonuclease: None detected.

Incubation of 54 ug (~ 32,000 units) with 110 ug of 43-mer oligonucleotide in 50 ul 50 mM Tris-HCl, pH 7.5, 5 mM MgCl<sub>2</sub>, 5 mM dithiothreitol, at 37 °C for 3 hours resulted in no observable degradation of the 43-mer oligonucleotide.

(3) Double Stranded Exonuclease: None detected.

Incubation of 54 ug (~ 32,000 units) with 50 ug of a 20 base pair double stranded DNA ladder in 50 ul 50 mM Tris-HCl, pH 7.5, 5 mM MgCl<sub>2</sub>, 5 mM dithiothreitol, at 37 °C for 3 hours resulted in no observable degradation of the 20 base pair dsDNA fragment, nor any of the upper DNA ladder bands.

(4) SDS-PAGE

SDS-PAGE analysis followed by coomassie blue staining shows only one major band. Purity is visually estimated at 98%. This product does not contain bovine serum albumin, or any other protein additives.

Made in USA by Edward Hyman

**BAYOU BIOLABS, LLC**

4724 Hessmer Avenue, Metairie, LA, 70002, USA  
phone: 504-723-1703 email: support@bayoubiolabs.com