

M13mp18 Single Stranded Phage DNA

Catalog # P-107

Lot: 21

Quantity Supplied: 500 ug

Concentration: 1.00 ug/ul

Storage Buffer: 10 mM Tris-Cl, 1 mM EDTA, pH 8.0

Store at -20 degrees C or below

Quality control:

0.5 ug phage DNA preparation was analyzed by agarose gel electrophoresis and stained by ethidium bromide. The gel demonstrated phage co-migration with control phage, an absence of RNA and very little chromosomal DNA contamination. Quantitation of phage DNA was performed by measuring the absorbance at 260 nm.

Recommendations:

Repeated freeze-thaw cycles have not been tested on this product. Vortex gently after thawing. This product may be diluted in TE buffer if it is too concentrated for routine use.

Made in USA by Edward Hyman

BAYOU BIOLABS, LLC

4724 Hessmer Avenue, Metairie, LA 70002, USA

phone 504-723-1703

support@bayoubiolabs.com

www.bayoubiolabs.com